

ABUS System Consumable Supplies

GE Healthcare offers consumable supplies for the ABUS (Automated Breast Ultrasound) systems including scan membranes and POLYSONIC® ultrasound lotion.

ABUS Single-Use Stabilization Membrane for the Invenia™ ABUS and somo•v™ ABUS Platinum

The ABUS disposable scan membrane provides comfort to the patient during the 3D image acquisition. The flexible, polyester, mesh fabric stabilizes the breast tissue and enables better image quality during the scan. The membrane is single-use-only which eliminates the risk of cross-contamination and the need for sterilization.

Curved membrane dimensions: 229mm (L) x 213mm (W) x 19mm (D)

| Catalog # | Description | Quantity per Box |
|-----------|---|------------------|
| E8340AB | ABUS Single-Use Stabilization Membrane for Invenia™ ABUS and somo•v™ ABUS Platinum models (non-Sterile) | 50/box |



PolySonic Lotion

POLYSONIC® Ultrasound Lotion is an acoustic coupling lotion that is hypo allergenic and has a light pleasing scent.

| Catalog # | Description |
|-----------|--|
| E8365EC | PolySonic Ultrasound Lotion (Non-Sterile) – 8.5 oz dispenser bottles (12 bottles/box) |
| E8365EF | PolySonic Ultrasound Lotion (Non-Sterile) – 1 US gallon with refillable dispenser (4 gallons/case) |



To reorder ABUS membranes:



Contact your local GE Healthcare representative.

To reorder PolySonic lotion:

Contact your local GE Healthcare representative.

Or

For a local distributor contact Parker Laboratories at:
Phone: (800) 631-8888
Email: parker@parkerlabs.com



To reorder Tristel Duo™ or Distel disinfecting wipes

Please contact your local Tristel representative.

Or

Tristel Headquarters at:
Phone: +44 (0)1638 721500
Email: mail@tristel.com



General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue any products described at any time without notice or obligation. Please contact your GE representative for the most current information.

©2016 General Electric Company.

GE and GE Monogram are trademarks of General Electric Company.

POLYSONIC is a trademark of Parker Laboratories, Inc. Tristel Duo is a trademark of Tristel Solutions Limited.

GE Healthcare, a division of General Electric Company.

GE Healthcare
3000 North Grandview
Waukesha, WI 53188
USA
www.gehealthcare.com



imagination at work